

Economy IV Pole



HCPCS:




FIGURE A

- 4 rubber wheels provide for a smooth transport and maneuverability
- Removable hooks make this economy IV pole versatile
- Easily converted to 2-hook or 4-hook with easy-to-release pushpin
- Chrome plated steel with weighted base provides strength, durability and reduces risk of tipping
- Height range of 40"82"
- Locking collar allows for easy height adjustment

PRODUCTS

Description

| | |
|---------|--------------------------------|
| 13029 | IV Pole,Econ.Rem.4Hook,1/cs |
| 13033 | IV Pole,Econ.Rem.2 hook,1/cs |
| 13033SV | IV Pole,Econ.SV,Rem.2Hook,1/cs |

SPECIFICATIONS

Description

| | |
|----------|------------------|
| Warranty | Limited Lifetime |
|----------|------------------|