drive DeVilbiss



Swivel Wheel with Lock & Two Sets of Rear Glides

HCPCS:



FIGURE A

FIGURE C

- Can be locked for fixed wheel or unlocked for a swivel wheel (Figure A B)
- Allows for 8 adjustments
- Rubber wheels allow walker to roll easily and smoothly over irregular surfaces
- Comes with rear glide caps, (item #10107) and glide covers (item #10107C), allowing use on all surfaces
- Weight Capacity 350 lbs

| PRODUCTS | Description |
|-----------------|----------------------------------|
| 10113 | Swivel-Lock Wheel,3",1pr/bx, RTL |
| 10115 | Swivel-Lock Wheel,5",1pr/bx, RTL |
| SPECIFICATIONS | Description |
| Weight Capacity | 350 lbs |