

Pulmo-Aide® Compact Compressor Nebulizer System

HCPCS: A7003,E0570



FIGURE A

- The Pulmo-Aide Compact Compressor Nebulizer System was built to last with high-quality metal mechanical parts, ensuring quality therapy for an extended period of time
- Small footprint
- Small, compact design
- Convenient carrying handle

PRODUCTS	Description
3655D	Pulmo-Aide Compact w/ Disp. Nebulizer
SPECIFICATIONS	Description
Warranty	Limited 5 year
Dimensions	7.5" (W) x 4.0" (H) x 7.2" (D)
Electrical Requirements	2.5 amps, 115 VAC 60 Hz
Free Air Flow	13 lpm Typical
Power Consumption	140 Watts Maximum
Weight	3.8 lbs
ACCESSORIES	Description
3655D-601	FILTERS (PKG. 5), EA
3655D-602	FILTERS (PKG. 25), EA
DL1050	PEDIATRIC CHARACTER MASK (50-PK), CS
3655D-621	VixOne Disp. Nebulizer, CS 50, CS