



HCPCS:







FIGURE A

FIGURE B

- High quality therapeutic foam mattress provides pressure redistribution and shear/friction reduction
- Drive's premium therapeutic Gravity foam mattress maximizes value and comfort with multiple layers of latex free foam
- Full top layer of open cell, visco-elastic memory foam improves immersion and pressure redistribution
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress
- High density foam perimeters provide solid edge support and facilitate safe patient transfers
- Available with or without raised side rails
- Heel section is sloped with channel cut foam center offloading pressure from vulnerable heels
- Removable multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS	Description
15970	Gravity 9 Mattress 80", 1/EA
15970-4280	Gravity 9 Mattress, 42"W x 80"L, 1 EA
15970-4284	Gravity 9 Mattress, 42"W x 84"L, 1 EA
15977	Gravity 9 MattressW/Elev Perim 80", 1/EA
15984	Gravity 9 Mattress 84", 1/EA
15996	Gravity 9 Mattress 76", 1/EA
G9-39X80X6	Gravity 9 Mattress 39"W x 80"L x 6"H
SPECIFICATIONS	Description
Warranty	36" W -15 Years Limited Foam 3 Year Cover 42" W- 7 Years Limited Foam 3 Year Cover
Side Rails	Yes
Cover Material	Stretch, Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
ACCESSORIES	Description
15977C	Mattress Cover 15977C,1/ea, EA
15970C	Mattress Cover 15970C,1/ea, EA
15984C	Mattress Cover 15984C,1/ea, EA
15996C	Mattress Cover 15996C,1/ea, EA