

## **Quick 'N Easy Comfort Mattress**

HCPCS: E0272-1,E0272









FIGURE A

FIGURE B

FIGURE C

- Durable high quality latex free polyurethane foam mattress for utmost comfort and support
- Non-slip bottom and 4 corner straps keep the mattress securely in place
- Two side carry handles make lifting and placement easier
- Compression technology creates a 14" x 14" x 37.5" box footprint creating convenience for transport and delivery
- Removable Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS	Description	
15076	Quick N Easy Comfort Mattress,1EA	
15076C	Cover,15076,1ea	
SPECIFICATIONS	Description	
Warranty	1 Year Limited	
Weight Capacity	275 lbs.	
	CFR 16 Part 1633 Federal Fire Code	
Cover Material	Fluid-Resistant Nylon	
Dimensions	36" (W) x 80" (L) x 6" (H)	
PARTS	Description	
15076C	Cover,15076,1ea, EA	