

Toilet Safety Rail

HCPCS:



FIGURE A



FIGURE B

- Allows individuals to safely sit or stand
- Easy, tool-free assembly and installation
- Comes with soft, padded handles for comfort
- Width adjusts from 16.5" - 19.5" between handles

PRODUCTS

Description

RTL12087	Toilet Support Rail,KD,RTL,1EA,2/CS
----------	-------------------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	24 lbs.
Warranty	Limited Lifetime
Weight	10.5 lbs.
Weight Capacity	300 lbs
Width Between Arms	16.5" - 19.5"