

## Steerable Knee Walker

HCPCS: E0118



- Knee walker can be steered for increased maneuverability
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot
- New, two-piece padded leg pad provides comfort
- Deluxe braking system
- Removable front basket
- Durable, steel frame with silver vein finish
- Tool free height adjustment
- 7" casters are ideal for indoor/outdoor use

**PRODUCTS****Description**

796	Steerable Knee Walker, Split KneePad, 1/EA
-----	--

**SPECIFICATIONS****Description**

Carton Shipping Weight	27 lbs
Handle Height	31"-40"
Knee Pad Dimensions	7" (W) x 15" (D)
Knee Pad Height	16.5"-20.5"
Unit Dimensions	16"-30"
Warranty	Limited Lifetime, Lifetime Warranty on Brake Cable
Weight	22 lbs
Weight Capacity	350 lbs

**PARTS****Description**

10230SN	Blow Molded Seat, 10230, 1ea, EA
9503H79014	Knee Pad, 790, 1ea, EA
95012R79625	Allen Bolts, 796, 4/pkg, PKG
9503H79021	Wheel Bearings, 790, PhoenixHD3, 1/pr, PR
95012R79601	Screws For Knee Pad Replacement, 796, 4/pkg, PKG
95012R79602	Brake Cable Only, 796, 1/ea, EA
95012R79603	Wheel Lock Without Pad, Left, 796, 1/ea, EA
95012R79604	Wheel Lock Without Pad, Right, 796, 1/e, EA
95012R79605	Wheel Lock Spring Only, 796, 1/ea, EA
95012R79606	Brake Pad Connector Bar W/Adj. Bolt, 796, 1/e, EA
95012R79607	Tiller Adjustment Knob, 796, 1/ea, EA
95012R79608	Axle Bolt, 796, 1/ea, EA
95012R79609	Front/Rear Wheel W/Bearings & Spacer, 796, 1e, EA

## PARTS

## Description

95012R79610	Handbrake W/ Cable, 796, 1/ea, EA
95012R79611	Knee Post Without Pad, 796, 1/ea, EA
95012R79612	Knee Pad Only W/Hardware, 796, 1/ea, EA
95012R79613	SeatPostReleaseLeverW/Hrdwr, 796, 1/ea, EA
95012R79614	Knee Post Adjustment Pin, 796, 1/ea, EA
95012R79615	Brake Pad W/Screws & Nuts, 796, 1/pr, PR
95012R79616	Fork Stem Bearing Assembly, 796, 1/pr, EA