



Bariatric Aluminum Folding Walker, Two Button

HCPCS:



FIGURE A

- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design accommodates individuals up to 600 lbs
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability
- · Vinyl-contoured hand grip

PRODUCTS	Description
10220-1WW	Walker, AdltHD, 2Btn w/Whls, 1/cs
SPECIFICATIONS	Description
Carton Shipping Weight	24 lbs.
Height	32" - 39"
Overall Width	27"
Warranty	Limited Lifetime
Weight	11 lbs.
Weight Capacity	600 lbs
Width	20.5" inside hand grip 23" inside back legs