

ENCORE[®] Sensi-Touch[®] PF

S U R G I C A L G L O V E S

Latex powder-free surgical gloves with a strong formulation, for everyday hand protection

- **Secure personal protection:** ENCORE[®] Sensi-Touch[®] PF surgical gloves are powder-free and feature a strong glove formulation, granting everyday protection you can depend on for general surgical procedures
- **Reliable comfort and fit:** These sterile surgical gloves have been designed using Ansell's proprietary anatomic formers to offer comfort and dexterous hand movement
- **Effortless donning:** Thanks to their polymer-coated internal finish, these dependable surgical gloves are both fast and easy to don
- **Reliable grip:** These latex surgical gloves also feature a smooth external surface, with a micro-textured finish, for secure instrument handling
- **Chemo-tested:** They have been tested for use with chemo-drugs, in accordance with ASTM D6978 requirements

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established - this product contains natural rubber latex, which may cause allergic reactions

Recommended For

- Brief, low-risk examination procedures
- Sterile procedures in the O.R.
- Use as an outer glove when double gloving
- Sterile procedures outside of the O.R.



ENCORE[®]
The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.

ENCORE[®] Sensi-Touch[®] PF

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Natural
Thickness	Standard
Cuff Style	Beaded
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	ASTM D3577, REACH, EN 420:2003 + A1:2009
Quality / Environmental Standards	ISO 13485, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	7821PF	7822PF	7823PF	7824PF	7825PF	7826PF	7827PF	7828PF

Contact your Ansell representative for ordering or more information.

North America
 US Tel: 800 952 9916
 CA Tel: 844 494 7854
 Email: insidesalesus@ansell.com
 Sandel US Tel: 866 764 3327
 Sandel email: sandel-cs@ansell.com

Central & South America
 Tel: +52(442) 296 20 50
 Email: cslac@ansell.com
Brazil
 Tel: +55-11-3356-3100
 Email: latam.medical@ansell.com

Europe, Middle East & Africa
 Tel: +32 (0) 2 528 74 00
 Email: info@ansell.eu
Asia Pacific
 Tel: +603 8310 6688
 Email: apac.medical@ansell.com

India
 Tel: +91 22 38124500
 Email: info.india@ansell.com
Australia & New Zealand
 Tel: +61 3 9270 7270
 Email: protection@ap.ansell.com
China
 Tel: +86(21) 3827 5005
 Email: infochina@ansell.com