



FreedomAlert[®]
Enjoy the Freedom of Independent Living

Active Landline Required (POTS, VoIP, DSL)

Water Resistant The pendant can be taken into the shower, where many serious accidents happen.

Voice Prompted System Provides announcements for set-up, battery test and full system check.

Pendant Battery A lithium-ion battery provides over 1 hour of talk time on a full charge, and 4 to 6 months in stand-by. Two rechargeable batteries are included. One charges in the base station charger, while the other is in the pendant.

System has standby emergency battery power Rechargeable batteries continually charge in the base station to provide up to 24 hours of standby power in the event of a power failure.

Multiple Pendants A total of up to 4 Pendants and Emergency Wall Communicators can be set to work with one system. Great for couples or assisted living facilities.

Pendant can be carried in 3 convenient ways. Included with each system is a **Lanyard**, a **Wrist Strap/Wheelchair holder** and a **Belt Clip**.

DECT Technology: Operation at 1.9 Ghz gives super clear voice quality as well as exceptional range.

Approved in all 50 states and Canada. FCC & IC certified. Factory is ISO-9001 approved. Conforms with UL1637.

LogicMark
The Leader in Two-Way Voice Communication


www.LogicMark.com

Pendant

Base Station

Optional Accessory:

With majority of falls occurring in the bathroom or bedroom, the **Emergency Wall Communicator** can be mounted anywhere to provide additional security when and where you need it.



1 Year Limited Warranty
FreedomAlert is protected
by US Patent #8,369,821



Brochure #6102
Copyright 2014 LogicMark, LLC V4.1

ENJOY THE FREEDOM OF INDEPENDENT LIVING



Family is just SECONDS away!

FreedomAlert[®]
Enjoy the Freedom of Independent Living

World's First Programmable
2-Way Voice Emergency
Pendant Communicator



No Activation Costs
No Contracts

Talk

Listen



LogicMark
The Leader in Two-Way Voice Communication



Complete coverage in and around your home

Because you are speaking through your pendant, FreedomAlert allows you to communicate your message from every corner of your home, outside in the yard and in the driveway. Emergencies can happen anywhere. And just one button activation.



FreedomAlert's 2-Way Voice Communication Range



Typical Personal Emergency Response System's 2-Way Voice Communication Range (speakerphone is in the base station)

The world's first 2-way emergency pendant communicator

Emergencies can happen in all parts of a home and frequently right outside a home in the yard or driveway. So whether you are just checking in with a friend, need some quick help from a neighbor, or are in an emergency, your FreedomAlert pendant will allow you to communicate your message instantly and efficiently from anywhere in and around your home. All with ONE button activation.

Up to 4 user programmable contacts to call anytime

You can personalize the system with a flick of a three position switch on the base unit to:

1. Call up to 4 custom contacts
2. Call up to 4 custom contacts & then the 911 emergency operator
3. Or just the 911 emergency operator.

It's your choice.

Emergency operators are available as backup if needed

FreedomAlert can also contact a 911 emergency operator as a backup in case you cannot reach a family member or neighbor. This way you are completely protected 24 hours a day, 7 days a week. You are never alone with FreedomAlert.

Traveling?

Whether you are moving across town or going to visit the grandkids across the country, this is not a problem for FreedomAlert. Just unplug your unit from the phone jack and toss it in your suitcase. When you get to your destination, simply plug it into the local phone jack. No additional programming needed. And if you call 911, they will have the local 'caller ID' of that location.

Never pay a monthly, activation or set up fee again.

FreedomAlert is yours to keep and use forever. Once you purchase it there is nothing more to pay.