



Thermoskin Hinged Elbow features two adjustable locking straps. Combining these straps with the medial and lateral hinges provides a bracing mechanism which stops hyperextension at zero degrees. This support may be suitable for those suffering from medial elbow instability or elbow hyperextension ailments.

PRODUCT ADVANTAGES:


- Lined with the exclusive Trioxon® Advantage material for heat & enhanced recovery
- Clinically tested to increase skin temperature and blood flow for pain relief
- Provides light but firm compression to counteract tissue swelling
- Adjustable locking straps provide additional support to prevent elbow hyperextension
- Universal design fits left or right elbow



THERMOSKIN®

TRIOXON® ADVANTAGE

Exclusive Trioxon 3 dimensional lining provides heat therapy & pain relief. Soft, comfortable inner layer wicks away moisture. Middle layer helps increase skin temperature for heat therapy. Outer layer provides light compression.



SIZING GUIDE

Measure around elbow joint with arm extended

PART NUMBER	SIZE	INCHES	CM
-	-	-	-
83137	S	9-10¼"	23-26
84137	M	10½-11¾"	27-30
85137	L	12-13¾"	31-35
86137	XL	14-15¾"	36-40
87137	XXL	16-17¾"	41-45