


CPAP Mask Accessories



Item #	Product Description	UOM
18217	Universal Chin Strap for CPAP	1/ea
C10200-014	Headgear for Full Face ComfortFit Deluxe EZ 100FDE Series	1/ea
C10200-018	Headgear for Full Face ComfortFit Deluxe 100FD Series	1/ea

 Deluxe, breathable, moisture wicking foam material, which keeps skin cool and comfortable increasing patient compliance

Universal Oxygen Cylinder Carry Bag for Wheelchair









- Adjustable, for use with "D" and "E" cylinders.
- Top and bottom clasps easily attach for secure fit to wheelchair.
- 2 carry handles.

Item #

STDS6008-1

- For use with most manufacturers' wheelchairs.
- Made of durable, easy-to-clean nylon.
- Size: 26.5" (H) x 7.5"(W) x 4.75" (D).

Universal Nebulizer Shoulder Carry Bag







Item #	Product Description	UOM	
18031	Universal Nebulizer Carry Shoulder Bag	1/ea	

- Designed to carry and store nebulizers.
- Universal size can accommodate all Drive nebulizers.
- Made of durable, lightweight, easy-to-clean nylon.
- Comes with an adjustable carry strap.
- Has a zippered cover and zippered pouch to conceal nebulizer accessories.
- Size: 9" (H) x 12" (W) x 7" (D).

CHAD® 3-in-1 Oxygen Cylinder Shoulder Carry Bag





Item #	Product Description	UOM
OP-150-800	CHAD® 3-in-1 Oxygen Cylinder Shoulder Bag	1/ea

- Hook and Loop Fastener brand touch storage pouch
- Opening for threading cannula on side of bag.
- Vertical, horizontal, or backpack configuration.
- Compatibility: M4, M6, M7, M9, while using any Drive or CHAD® conservers and regulators.
- Size: 15.75" (H) x 4.33" (W) (at base) x 8.00" (W) (at widest point).