


## The IontoPatch Family of Products

Four products provide different dosage and treatment options based upon each patient's clinical presentation and tolerance.




**NEW!**

### IontoPatch **EXTRA STRENGTH**®

Stronger dose (120mA-min) for more aggressive yet safe treatment. Treats most areas, including elbows, knees, wrists, feet and shoulders. 8-hour average patient wear time.\*


- ▶ More therapy in fewer treatments.
- ▶ Patients with severe conditions improve faster allowing a full treatment regimen sooner.



### IontoPatch® 80

Treats most areas, including elbows, knees, wrists, feet and shoulders. 14-hour average patient wear time.\* 80mA-min dose strength


- ▶ Lowest current to treat patients with sensitive or fair skin, including the elderly.



### IontoPatch **STAT**®

The same IontoPatch 80 benefits in shorter 4-hour average patient wear time.\* Treats most areas, including feet, elbows, knees, wrists, and shoulders. 80mA-min dose strength

- ▶ Treats patients who require or desire shorter wear time.



### IontoPatch® SP

Treats smaller areas, including fingers and Achilles tendon. 14-hour average patient wear time.\* 40mA-min dose strength

- ▶ Designed to adhere and contour to small anatomical areas like fingers and Achilles.

\*Delivery time is dependent upon individual patient skin resistance. Study data on file.

- Appropriate for both acute and chronic conditions
- Family of four products that treat large and small areas, offering three different wear times and three different dose strengths
- Wrappable, compressible, weight bearing
- The Original Iontophoretic Patch proven by experience:
  - over 6 million administered.
  - supported by peer-reviewed published studies.
  - first and only independent published outcome-based study comparing iontophoresis to injection.
- No charging station or controller required
- Extended time-released iontophoresis
- More cost-effective and less time consuming
- More drug delivered and remains localized



**Patterson Medical product code:**  
IontoPatch Extra Strength® 5509-36  
IontoPatch® 80 A939-65-0  
IontoPatchSTAT® 5640-91  
IontoPatch® SP 9226-40



See the link to the IontoPatch Application Instruction Video located at [www.iontopatch.com](http://www.iontopatch.com) or on the IontoPatch page of the Patterson Medical website at [www.PattersonMedical.com](http://www.PattersonMedical.com)