

MAKES ORDINARY DONUT RINGS OBSOLETE

3-in-1 Complete Comfort Seating Solution

The Kabooti® Seat Cushion makes other donut pillows obsolete. It's a 3-in-1 complete comfort seating solution – coccyx cushion, seating wedge and donut ring all in one – to provide the perfect posture, pressure-free comfort, and prolonged seating relief.

- Generous Coccyx Cutout
- Contoured for Comfort
- Tapered Front Edge
- Reduces Lower Back Pain
- Provides Stabilization



Guard your Kabooti® Seat Cushion against incontinence, perspiration, spills and more with the Kabooti® Waterproof Cover.

(fits Kabooti® Original and Kabooti® Ice)

Sold Separately

Transform your Kabooti® into a Kabooti® Ice with this Icy Gel Insert. This reusable ice pack provides cooling relief for sensitive areas.

Sold Separately



Fun Name! Award-Winning Package Design Merchandises Beautifully!




Also available in Kabooti® Large - 18% Wider!


KABOOTI®

THE BETTER DONUT

Ergonomic Design



3-in-1 Complete Comfort Seating Solution



Coccyx Cushion




Seating Wedge



Donut Ring

The KABOOTI® is uniquely engineered to be superior to conventional donut cushions.



Pressure mapping* shows a comfortable seating with minimal pressure on the Kabooti, but the donut cushion shows extremely high pressure concentration (red and orange areas) on buttocks area. This can result in uncomfortable seating, reduced skin and surface blood circulation and potential damage to skin.

*Pressure maps performed with 175 lb. male, using XSensor pressure mapping equipment & software.

Product	Item #	UPC	Prod Dims (in) L x W x H	Pckg Dims (in) Lx W x H	Pckg Wt. (lbs.)	Case Pack	Master Ctn (in) L x W x H	Carton Wt. (lbs)	Cartons/ Pallet
Kabooti® - Blue	30-751RB	737709006311	17.5 x 13.5 x 3	17.5 x 16 x 3	1.6	4	17.5 x 14.5 x 10.3	7.2	48
Kabooti® - Gray	30-751RGY	737709003628	17.5 x 13.5 x 3	17.5 x 16 x 3	1.6	4	17.5 x 14.5 x 10.3	7.2	48
Kabooti® Large - Blue	30-740RB	737709004472	20 x 14 x 3	20 x 15.5 x 3	2	4	21 x 14.5 x 10.5	7.8	48
Kabooti® Large - Gray	30-740RGY	737709004489	20 x 14 x 3	20 x 15.5 x 3	2	4	21 x 14.5 x 10.5	7.8	48
Kabooti®Ice	30-761RBI	737709004816	13.5 x 17 x 3	15.5 x 17 x 3	2.2	4	17.5 x 14.5 x 10.3	8.7	48
Kabooti® Icy Gel Insert	30-761V2GEL	737709005417	6.5 x 5 x 2	7 x 4 x 7	1.5	4	15 x 21.5 x 13	6.2	
Kabooti® Waterproof Cover	1-634WP-100R	737709005103	17.5 x 14 x 3.5	6.5 x 1 x 10	2.9	4	7 x 9.5 x 4.5	1.4	



DME/HME Division
855-240-3301
sales@contourhealth.com

www.contourproducts.com

Retail Division
800-692-6686
sales@contourhealth.com



Invented and designed in the United States. Made in China. Contour is a registered Trademark used by Contour Products, Inc., under license from Banyan Licensing, L.L.C. © 2020. All rights reserved.

CO037_1220