

JEANIE RUB[®] VARIABLE SPEED BODY MASSAGER

- Deep, penetrating massage with a twist of the switch
- Helps reduce tension and rejuvenate sore, tired muscles
- Promotes circulation to increase blood flow
- Helps increase elasticity of muscles
- Trusted by healthcare professionals
 worldwide
- Variable speeds provide soothing massage or more invigorating massage
- Provides the luxury of a full body massage without the recurring cost of massage therapists



OTHER SPECIAL FEATURES

- Built to ANSI standards; CSA certified
- 500-4600 RPM orbital massaging action
 10 feat neuron const
- 12-foot power cord
- Requires 120 VAC Nominal Power Source

Available From:



LUXURIOUS MASSAGE

The Jeanie Rub Massager[®] with variable speed provides a deep, penetrating massage with the twist of a switch. Help reduce pain, muscle tension, and stiffness by adjusting massage intensity from a soft, gentle speed to a deep, invigorating speed. Treat yourself to the luxury of a full body massage without the recurring cost of a massage therapist.

The soothing, deep penetrating massage increases blood flow to assist in rejuvenation of sore, tired, or tense muscles. Adding the Jeanie Rub Massager to your pre-workout routine can help increase elasticity of the muscles to help prevent workout related injuries.

VARIABLE SPEEDS

This massager delivers an invigorating 500-4600 RPM orbital massaging action; lower speeds deliver a soothing massage while faster speeds are more invigorating. The orbital action of the Jeanie Rub provides a softer and deeper massage than vibrating or percussion type massagers. This product is built with a cast aluminum construction with a black powder coat finish that looks great and resists abrasions, ensuring years of use.

SPECIFICATIONS

Product Selection: Jeanie Rub® Variable Speed Body Massager, Black PRO-3401

Product Weight: 7 lbs

Dimensions: 10.25" x 8.25" x 4.625" Pad measures: 4" x 9"

Warranty: 1 Year

Case Quantity: 4

Case Dimensions: 19.1215" x 10.875" x 8.5"

Case Weight: 28 lbs

Latex: This product is not made with natural rubber latex

Country of Origin: Assembled in the USA of globally sourced materials

California Residents Only: **WARNING:** Cancer and Reproductive Harm – www.P65Warnings.ca.gov



JEANIE RUB[®] MASSAGER ACCESSORIES

All massager accessories are:





Wall Hanger ACC-884

Product Dimensions: 14" x 5" x 5" Product Weight: 1.5 lbs Materials: Mixed Hardwoods Case Quantity: 1 Case Weight: 3 lbs Case Dimensions: 10.25" x 13.5" x 12.5"



PROFESSIONAL KIT INCLUDES

- Jeanie Rub Massager | PRO-3401
- Paraspinal Accessory | ACC-896
- Extremity Accessory ACC-897
- Carrying Case | ACC-885-BK-JRM
- Fleece Pad Cover | ACC-882

Available From:



Fleece Pad Cover ACC-882

Product Dimensions: 9.5" x 5.5" x .5" Product Weight: .1 lbs Materials: 100% polyester Case Quantity: 12 Case Weight: 4 lbs Case Dimensions: 11.25" x 9" x 6.25"



Sheepskin Pad Cover ACC-886 Product Dimensions: 10" x 8" x 1" Product Weight: 3 lbs Materials: Australian Sheepskin Case Quantity: 6 Case Weight: 3 lbs Case Dimensions: 11.25" x 9" x 6.25"



Sure-Chek[®] Covers - 3 Pack ACC-892-GR

Product Dimensions: 9" x 4" Product Weight: .15 lbs Materials: Sure-Chek, elastic Case Quantity: 3 Case Weight: 2 lbs Case Dimensions: 11.25" x 9" x 6.25"



Extremity Accessory ACC-897

Product Dimensions: 10.25" x 5" x 4.5" Product Weight: .2 lbs Materials: Polyethylene foam Case Quantity: 6 Case Weight: 2 lbs Case Dimensions: 11" x 14" x 13"

©2022 Core Products International, Inc. / Osceola, Wisconsin



Nylon Carrying Case ACC-885-BK-JRM Product Dimensions: 12.5" x 8" x 12.5" Product Weight: 1.05 lbs Materials: Nylon, Foam, Polyester, pvc coating Case Quantity: 6 Case Weight: 8 lbs Case Dimensions: 16" x 14.25" x 13.5"



Multi-Use Holder ACC-887 Product Dimensions: 13" × 8.5" × 2.25" Product Weight: 4 lbs Materials: Polyurethane Foam, Volextera Case Quantity: 5 Case Weight: 4 lbs Case Dimensions: 11" × 14" × 23"



Terry Cloth Covers - 3 Pack ACC-893 Product Dimensions: 10" × 7"

Product Dimensions: 10 x 7 Product Weight: .2 lbs Materials: 10oz knit terry: 70% Cotton/30% Polyester Case Quantity: 3 Case Weight: 2 lbs Case Dimensions: 11.25" x 9" x 6.25"



Paraspinal Accessory ACC-896 Product Dimensions: 10.25" x 5" x 2" Product Weight: .15 lbs Materials: Polyethylene foam Case Quantity: 6 Case Weight: 2 lbs Case Dimensions: 10.25" x 13.5" x 12.5"