









# Nebulizer Compatibility - Which Nebulizers work with which Compressors?






Customers often have a question as to which nebulizer will work with their compressor. Can they use a different brand's nebulizer with their compressor?

So, the answer is that the nebulizer sets are interchangeable. Reusable and disposable nebulizers are not manufacturer or brand specific, but are all compatible with every brand. They only require compressed air from the compressor machine to work. As long as the tubing fits both the nebulizer and compressor, your machine is set to work. And to make this possible, usually the tubing size is made universal. But still it is the best to read product details to confirm the compatibility.

Product Image	Product Name	Product Type
	<p><b>Drive Disposable Nebulizer Kit</b></p> <p>#Item - NEBKIT500, NEBKIT600</p> <p><a href="#">View Details</a></p>	Disposable
	<p><b>Drive Reusable Nebulizer Kit</b></p> <p>#Item - 18254</p> <p><a href="#">View Details</a></p>	Reusable
	<p><b>Drive Medel JetNeb Plus Reusable Nebulizer</b></p> <p>#Item - 18260</p> <p><a href="#">View Details</a></p>	Reusable

	<p><b>Devilbiss VixOne Disposable Nebulizer</b></p> <p>#Item - 3655D-621</p> <p><a href="#">View Details</a></p>	<p>Disposable</p>
	<p><b>Respironics Sidestream Custom Disposable Nebulizer Kit</b></p> <p>#Item - 1028125</p> <p><a href="#">View Details</a></p>	<p>Disposable</p>
	<p><b>Respironics SideStream Disposable Nebulizer</b></p> <p>#Item - HS800</p> <p><a href="#">View Details</a></p>	<p>Disposable</p>
	<p><b>Respironics Sidestream Reusable Nebulizer</b></p> <p>#Item - HS860</p> <p><a href="#">View Details</a></p>	<p>Reusable</p>
	<p><b>Respironics SideStream Plus Breath Enhanced Reusable Nebulizer</b></p> <p>#Item - HS870</p> <p><a href="#">View Details</a></p>	<p>Reusable</p>

	<p><b>Pari LC Plus Reusable Nebulizer</b></p> <p>#Item - 022F81</p> <p><a href="#">View Details</a></p>	<p>Reusable</p>
	<p><b>Pari LC Star Reusable Nebulizer</b></p> <p>#Item - 022F51</p> <p><a href="#">View Details</a></p>	<p>Reusable</p>
	<p><b>Pari LC Sprint Reusable Nebulizer</b></p> <p>#Item - 23F35</p> <p><a href="#">View Details</a></p>	<p>Reusable</p>
	<p><b>CareFusion AirLife Misty Max 10 Disposable Nebulizer</b></p> <p>#Item - 2430, 2432, 2433, 2434, 2435, 2438, 2439, 2440, 2443, 2444, 2445, 2446</p> <p><a href="#">View Details</a></p>	<p>Disposable</p>
	<p><b>Carefusion Airlife Sidestream High-Efficiency Nebulizer</b></p> <p>#Item - 2173</p> <p><a href="#">View Details</a></p>	<p>Disposable</p>

	<p><b>CareFusion AirLife Infant Nebulizer Kit with Custom Misty Nebulizer</b></p> <p>#Item - 6200-504</p> <p><a href="#">View Details</a></p>	<p>Reusable</p>
	<p><b>Salter 8900 Series Disposable Small Volume Jet Nebulizer With Supply Tube</b></p> <p>#Item - 8900-7, 8906-7, 8911-7, 8924-7, 8900-7-50, 8906-7-50, 8911-7-50, 8924-7-50</p> <p><a href="#">View Details</a></p>	<p>Disposable</p>
	<p><b>Sunset Healthcare Disposable Nebulizer Kit with Mask</b></p> <p>#Item - RES092</p> <p><a href="#">View Details</a></p>	<p>Disposable</p>
	<p><b>Sunset Healthcare Reusable Nebulizer Kit with T-Piece</b></p> <p>#Item - RES093</p> <p><a href="#">View Details</a></p>	<p>Reusable</p>
	<p><b>Hudson RCI Micro Mist Small Volume Nebulizer</b></p> <p>#Item - 1882, 1883</p> <p><a href="#">View Details</a></p>	<p>Disposable</p>

	<p><b>Omron A.I.R.S. Disposable Nebulizer Kit</b></p> <p>#Item - 9911</p> <p><a href="#">View Details</a></p>	<p>Disposable</p>
	<p><b>Graham Field Complete Nebulizer Set</b></p> <p>#Item - GF61399</p> <p><a href="#">View Details</a></p>	<p>Reusable</p>
	<p><b>Allied Healthcare Inc Hand Held Nebulizer with Mouthpiece and Tee</b></p> <p>#Item - 61399, 61400</p> <p><a href="#">View Details</a></p>	<p>Reusable</p>
	<p><b>McKesson Nebulizer Kit Infant Pacifier</b></p> <p>#Item - 32641</p> <p><a href="#">View Details</a></p>	<p>Reusable</p>