

I N S T R U C T I O N S

Pistol-Grip Remote Toenail Clipper

Catalog Number 741390000

U.S. Patent #4,893,406

The Remote Toenail Clipper gives you an extra 18½" (47 cm) of reach that makes it possible to clip toenails without bending.

Before using, observe the operation of the clipper by squeezing the trigger and noting the motion of the clipping head.

To clip toenails, slip the clipping blades over the toenail and squeeze the trigger. Be careful not to pinch any part of the toe or to clip too deeply. Release trigger and continue on.

Keep out of reach of children.



I N S T R U C T I O N S

Pistol-Grip Remote Toenail Clipper

Catalog Number 741390000


U.S. Patent #4,893,406

The Remote Toenail Clipper gives you an extra 18½" (47 cm) of reach that makes it possible to clip toenails without bending.

Before using, observe the operation of the clipper by squeezing the trigger and noting the motion of the clipping head.

To clip toenails, slip the clipping blades over the toenail and squeeze the trigger. Be careful not to pinch any part of the toe or to clip too deeply. Release trigger and continue on.

Keep out of reach of children.



A Division of SP Industries

661 Route 23 South • Wayne, NJ 07470

973-628-7600 • email custservice@maddak.com • www.spableware.com

974139006

0716



A Division of SP Industries

661 Route 23 South • Wayne, NJ 07470

973-628-7600 • email custservice@maddak.com • www.spableware.com

974139006

0716