

# Crosstex™ Plastic Headrest Covers

**REF** LOCP / LOCPW / L3CP / L3CPW / L14CP

## INSTRUCTIONS FOR USE

- 1. Use a surface disinfectant as recommended by the chair manufacturer to clean the headrest.
- 2. Hold cover with opening on the bottom.
- 3. Place cover over headrest until headrest is completely covered.



**CROSSTEX**  
 Crosstex International, Inc.  
 6789 W. Henrietta Road, Rush, NY 14543 USA  
 585.359.0130 | 631.582.6777 | crosstex.com

**EC REP** MDSS GmbH  
 Schiffgraben 41  
 30175 Hannover, Germany **CE**

**MD**

Not made with natural rubber latex.

Made in USA  
 PLASTICHEADCVR-IFU-EN / REV C - 08/2021  
 All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.  
 ©2021 Crosstex International, Inc. All rights reserved.