



PROCARE®

WeekENDER™ Recreational Activity Brace

Lightweight, rigid, low profile. May be used for post-op, rehabilitation or functional application. Dual axis motion control hinges with foam/chamois condyle pads provide full and limited range-of-motion with flexion/extension stops at 0°, 15°, 30°, 60° and 90°.

Indication: Recommended for mild medial/lateral instability, including ACL/PCL laxity.

Right	Left	Circum.*	Size
79-94353	79-94363	14.5" – 18.5"	S
79-94355	79-94365	18.5" – 21"	M
79-94357	79-94367	21" – 23.5"	L
79-94358	79-94368	23.5" – 26.5"	XL
79-94359	79-94369	26.5" – 29"	2XL
Replacement Strap with Standard Foam Pads			
79-94373	Thigh, Calf & Condyle Pads		S
79-94375	Thigh, Calf & Condyle Pads		M
79-94377	Thigh, Calf & Condyle Pads		L
79-94378	Thigh, Calf & Condyle Pads		XL
79-94379	Thigh, Calf & Condyle Pads		2XL
Chamois Replacement Pad Set			
79-94323	Thigh, Calf & Condyle Pads		S
79-94325	Thigh, Calf & Condyle Pads		M
79-94327	Thigh, Calf & Condyle Pads		L
79-94328	Thigh, Calf & Condyle Pads		XL
79-94329	Thigh, Calf & Condyle Pads		2XL

*Measurement taken 6" above mid-patella.