



TECHNICAL DATA SHEET



PRODUCT Crosstex Isofluid®, Isofluid Fog-Free®, Isofluid Fog Free with Shield, Isofluid SecureFit®, Isofluid Fog-Free SecureFit, Isofluid Fog-Free SecureFit with Shield – Face Masks

- MATERIALS**
- Outer Material – 100% Spunbond polypropylene
 - Inner Material – Medical Grade Tissue
 - Filter material – 100% Meltblown polypropylene, 18 gsm
 - Earloop Material – Flat, white knitted Elastic, latex free
 - Nose/Chin piece – Aluminum strip
 - Binding material – 9/16" 1.25 oz spunbond polypropylene
 - Vapor Barrier – Fog Free masks ONLY, Polyethylene coated polypropylene 1.25", 1.0 oz
 - Anti Fog/Fog Free Strip – 100% meltblown polypropylene 1 1/6", 1.9oz
 - Shield – 4 mil mylar polyester film
 - Colors – Blue, Green, Lavender, Pink, Sapphire, Teal, Topaz, Turquoise & Kaleidoscope.

- DETAILS**
- Meets ASTM F2100-11 Standards, ASTM Level 1
 - Meets EN 14683 Type II Requirements

PRODUCT/ PACKAGING REQUIREMENTS	Isofluid – 50/Box	Item # GCIBL – Blue	Item # GCIGR – Green	Item # GCIPK – Pink
		Item # GCITE – Teal	Item # GCISA – Sapphire	Item # GCIBL – Topaz
		Item # GCILV – Lavender	Item # GCITQ – Turquoise	
	Isofluid Fog-Free – 40/Box	Item # GCICXB – Blue	Item # GCICXS – Sapphire	
		Item # GCICXZ – Topaz	Item # GCICXT – Turquoise	
	Isofluid Fog Free with Shield – 25/Box	Item # GCIPWB – Blue	Item # GCIPWS – Sapphire	
		Item # GCIPWZ – Topaz	Item # GCIPWT – Turquoise	
	Isofluid Secure Fit – 50/Box	Item # GCIBLSF – Blue		
		Item # GCILVSF – Lavender		
		Item # GCIPKSF – Pink		
	Isofluid Bulk Pack– 100/Box	Item # GCIBL100 – Blue		
	Isofluid Fog-Free Secure Fit	Item # GCICXBSF – Blue – 40/Box		
	Isofluid Fog-Free Secure Fit with Shield	Item # GCIPWBSF – Blue – 25/Box		

STORAGE REQUIREMENTS N/A

SHELF LIFE N/A

Crosstex®, Isofluid®, Isofluid Fog-Free®, and SecureFit® are registered trademarks of Crosstex International.