

# EVAC+SLED

THE EVAC+SLED ENABLES CARERS TO RAPIDLY AND SAFELY EVACUATE PEOPLE WITH MOBILITY PROBLEMS IN AN EMERGENCY.

The **Evac+Sled** is an evacuation board designed to help evacuate non-ambulant people or people who may not be able to walk down stairs in an emergency. It can be used for children to adults 6 foot 6 inches (198cm.) tall.

The **Evac+Sled** can be used in office buildings, theatres, libraries, hospitals, schools and many other facilities.



## Features

### 1 Long pulling straps

Located at both head and feet to allow for multiple users

### 2 Angled side board flaps

Insures the body is kept in place and offering protection when transported around doorways etc.

### 3 Secured patient safety straps

For increased passenger safety andn comfort

### 4 Adjustable pillow

For added head protection and allowing to fit persons of varying heights

### 5 Glide-able plastic board

Corrugated plastic board allows for excellent heat insulation preventing injury, plus removing the necessity to unroll every few months to prevent the plastic keeping its shape

## Includes

- 440lbs Weight capacity
- 2 person operation
- Can be moved easily over virtually any surface
- Narrow width and flexibility
- Wipe clean storage bag
- Can be used down stairs

## Benefits

- No installation necessary, ready to use instantly
- No belay system required
- No dangerous ropes to negotiate during evacuation



### Dimensions packed

Height	22"
Width	27"
Depth	5"
Weight	4.8lbs

### Dimensions unpacked

Height	89"
Width	51"
Depth	3"
Weight Capacity	440lbs

For over 40 years EVAC+CHAIR® is America's Original and World's Number One Selling Evacuation Chair. Don't Risk Lives, Beware of Poor Imitations and Inferior Substitutions. There is only one EVAC+CHAIR®.

Unsurpassed Quality and Workmanship, Successful in Saving Lives in 9/11, 2 YEAR MANUFACTURERS WARRANTY.

For further information please contact:

Evac+Chair North America LLC

3000 Marcus Ave, Suite #3E6, Lake Success, NY, 11042-1012

T 516-502-4240 W www.evac-chair.com E sales@evac-chair.com

# EVAC+CHAIR®

THE WORLD'S NO.1 STAIRWAY EVACUATION CHAIR SINCE 1982