

# BioMed® Revived II

Digital TENS / EMS / Massage Stimulator



"Live and Enjoy Life Pain Free"



# BioMed® Revived II



## PACKAGE INCLUDES

- 1 BioMed® Revived II
- 4 2"x 2" electrodes
- 2 Lead wires
- 3 AAA alkaline batteries
- 1 Instruction manual

## SPECIAL FEATURES

- Auto shut off & lock function
- Low battery status indicator
- No-load detection
- Programed & adjustable modes
- HCPCS if billed 0730, 0745

FDA REGISTERED CLASS II MEDICAL DEVICE



## GENERAL SPECIFICATIONS

SIZE	4.29" x 2.15" x .90"
WEIGHT	2.47 oz.
POWER SOURCE	3 AAA 1.5V alkaline batteries (do not use Lithium)
CHANNELS	Dual (2)
WAVEFORM	Symmetrical, biphasic square pulse Asymmetrical, biphasic square pulse
PULSE RATE	2 - 100 Hz (hertz) adjustable
PULSE WIDTH	100 - 300 µS (microseconds) adjustable
AMPLITUDE	0- 120 mA peak
INTENSITY	40 levels
TIMER	5 - 90 minutes
TOLERANCES:	+/- 10%

Data recorded across a 500 OHM load resistance

BioMedical Life Systems has been producing the best portable therapy devices and accessories for over thirty-five years. The result is a comprehensive family of products meeting the needs of physical therapists, pain clinics, medical insurance, veterinarians, professional athletics, and individuals the world over.

Does pain management have to involve drugs? Absolutely not! Pain can affect your mood, social events, and joyful activities with family & friends. Researchers agree that the use of an electrotherapy device can provide sufficient relief from neuropathy, muscle, or joint pain so you can live and enjoy life!

The BioMed® Revived II is a digital, lightweight, and fully functional dual-channel therapy stimulator.

- 6 preset and 3 adjustable TENS programs for pain.
- 5 preset and 3 adjustable EMS programs for muscle.
- 5 massage programs to relax. Incorporates a timer.

With 40 levels of intensity, a timer, and a manual for detailed program settings, your sure to find the perfect setting for your personnel therapy sessions.



1954 Kellogg Avenue  
Carlsbad, CA 92008  
Ph: (800) 726-8367 USA  
Ph: (760) 579-0801 INTL  
Fx: (760) 929-9953  
www.bmls.com  
information@bmls.com



Cage Code 3R8H6

