



Completely Coordinated.
Completely Cost Effective.
Completely Easy.

ZOO BUS SCALE COMPLETE



Fun Series, Zoo Bus, Pediatric Scale Table



MODEL 7822

TABLE FEATURES

- Clinton original illustrations
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Built-in concealed paper dispenser and cutter
- Adjustable built-in levelers provide a stable surface
- Easy-clean all laminate base
- Optional storage area locks available
- Ample 2-door storage
- Soft-close, door hinges

SCALE FEATURES

- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.



MODEL 7822 SPECIFICATIONS

Top:	43"L x 24"W
Base:	47"L x 21"W x 36"H

California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Upholstery



Gunmetal

Schoolhouse Cabinets

MODEL 6120-BW

"Schoolhouse" Base Cabinets

4 doors with adjustable shelves and 3 drawers.
Shipped in one unit.

"Schoolhouse" Wall Cabinets

4 doors with adjustable shelves. Shipped in one unit.

Base & wall cabinets are only sold in sets

- Clinton original illustrations
- Pre-assembled units save time
- Easy-clean all laminate construction
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Tough, chip-resistant PVC edges
- Assembled with steel fasteners
- 4 adjustable leveling feet on each unit
- 3-way adjustable, soft-close 110° hinges
- Flat panel drawer fronts with metal sides that hold up to 75 lbs.
- Recessed toe kick
- Easy hook and rail installation for upper cabinets
- Slate Gray edgebanding on all edges
- Customizable with 4 options
- Easily moves for relocation
- Drawers can be added to any base unit
- Doors and drawer can be moved to accommodate sink location
- Shown with optional sink and faucet



SPECIFICATIONS

BASE: 60"L x 18"D x 35"H

WALL: 60"L x 12"D x 24"H

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.


1



Options



066 Door & Inside Latch Combo



055 Door Locks




044 Add a Cabinet Drawer




22 Sink & Faucet

Optional Items in Red



1 Square = 1 foot

NOTE: These are suggested floor plans. Please check local building and fire codes before beginning any construction or renovations.



Paint option with chair rail



Walls

2022-40 Banana Yellow

795 Highland Blue

Chair Rail/Trim

797 Athens Blue

Flooring

Medintone® H5301 Gray Light

Medintone® H5302 Deep Gray

OPTIONAL ACCESSORIES *(Shown in photo on front)*



- 01-CC Koala Pyramid Wall Sticker (34" x 73")
- C-40 Side Chairs in Royal Blue upholstery (3RB)
- T-70 Clear Labeled Sundry Jars (5 in set)
- GS-3022 Double Slimline Stainless Steel Glove Box Holder



800-441-9131 • 717-848-3519

Toll-Free Fax 866-522-1233

Fax 717-843-5871

www.clinton-ind.com