


## Flat Top, Imaging Table with Drop Window



### MODEL 80071

- One piece flat top with drop window
- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser, cutter and positioning straps included
- Optional side rails & casters available
- Model 80071: 500 lbs load capacity under normal use
- Model 80071-X: 600 lbs load capacity under normal use

California residents please be advised of the following, pursuant to Proposition 65:



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



*Trendelenburg*



*Reverse Trendelenburg*



*Optional Side Rails & Casters*

### 80071 SPECIFICATIONS

	length	width	height range
overall	72"	27"	21 <sup>1</sup> / <sub>2</sub> " – 37 <sup>1</sup> / <sub>2</sub> "
window	10"	10"	

**80071-X** Extra-Wide 34" width model



### Upholstery Colors



## OPTIONS

### Side Rails (pair)

- 3 height positions
- Plunger height adjustment

098



### Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

