

Fashion Finish, Molded Top Treatment Cabinet with 2 Doors and 2 Drawers



tashiontinish

MODEL 8822-AF


- Cabinet with 2 doors, 2 drawer and 1 adjustable shelf
- Easy-clean, all-laminate interior
- Tough, Fashion Finish thermofoil exterior with no front-facing seams
- Laminate or optional Fashion Finish exterior
- One-piece, vacuum-formed, gray-haircell, textured, top rimmed on all sides
- Recessed toe-kick
- Fully adjustable, soft-close, door hinges
- Smooth glide, drawers with one-piece , metal sides and a load capacity of 75
- Optional locks are available

8822-AF SPECIFICATIONS

24.875 W x 18.375 D x 35" H



California residents please be advised of the following, pursuant to Proposition 65:



Options







066 Door and Inside Latch Combo

Finishes







Twilight

Arctic White M

Metropolis Gray