

Fashion Finish, Molded Top, Mobile Treatment Cabinet with 2 Doors and 2 Drawers



fashionfinish
CASEWORK BY CLINTON

MODEL 8922-AF


- Mobile cabinet with 2 doors, 2 drawers and 1 adjustable shelf
- Easy-clean, all-laminate interior
- Tough, Fashion Finish thermofoil exterior with no front-facing seams
- One-piece, vacuum-formed, gray-haircell, textured, top rimmed on all sides
- Fully adjustable, soft-close, door hinges
- Smooth glide, drawers with one-piece metal sides and a load capacity of 75 lbs.
- Tough, plastic bumpers on each side of the base
- Cabinet is mounted on four, 3" swivel casters (two locking)
- Optional locks are available

8922-AF SPECIFICATIONS

24.875 W x 18.375 D x 35.25" H



California residents please be advised of the following, pursuant to Proposition 65:




WARNING: This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: www.P65Warning.ca.gov.

Options



055 Door & Drawer Locks



066 Door and Inside Latch Combo

Finishes



Twilight



Arctic White



Metropolis Gray