


## \*Six-Leg, 4-Hook Stainless Steel Infusion Pump Stand



### MODEL: IVS-734

- 4-hook top
- Stainless steel pole/knob-lock adjustment
- 6-leg, powder-coated, metal base
- Oversized steel base for stability
- Detachable ram's horn hooks with plastic tips
- 4" heavy duty, rubber wheel, ball bearing casters
- Extra rigid, pole-in-tube mounting system
- Heavy duty, low center of gravity base
- See available options
- 65 lbs. load capacity under normal use

*\*Clinton stainless steel items are not M.R.I. compatible.*


**Shipping Weight: 18 lbs.**

#### SPECIFICATIONS

Base Diameter—23"

Height Range: 60"—100"

**California residents please be advised of the following, pursuant to Proposition 65:**

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

**OPTIONS:**



*IV-46 Pump Tray Support*



*IV-41 Universal Clamp*



*IV-55 Extra Hooks and Organizer Hooks*



*Heavy Duty Black Epoxy Wire Basket*  
*IV-51B - 14" x 8" x 8.5"*  
*IV-52B - 12" x 6" x 6.5"*



*Heavy Duty Stainless Steel Wire Basket*  
*IV-51S - 14" x 8" x 8.5"*  
*IV-52S - 12" x 6" x 6.5"*



*IV-57 Single Ram's Horn*



*IV-44 Oxygen Cylinder Holder*



*IV-43 Steering Handle*



*IV-56 Drainage Hook*