

Hurricane® Freedom Edition

HCPCS:



- Freestanding design
- Folds & unfolds in seconds
- 8 convenient height adjustments
- Whisperflex pivoting base design for stability & balance
- Steadigrip base for all-terrain traction
- Comfort-fit shock-absorbing handle
- Anodized aluminum shaft
- Wrist strap included
- Floor Display Contains: 4 Black, 2 Red, 2 Blue, 4 Purple

PRODUCTS

Description

HCANE-BK-C2	Hurricane,freedom,BLK,1/EA,6/CV
HCANE-BL-C2	Hurricane,freedom, BLUE,1/EA,6/CV
HCANE-PB-C2	Hurricane, freedom, Pastel Blue
HCANE-PG-C2	Hurricane, freedom, Pastel Green
HCANE-PP-C2	Hurricane, freedom, Pastel Pink
HCANE-PR-C2	Hurricane,freedom, PURPLE,1/EA,6/CV
HCANE-RD-C2	Hurricane,freedom, RED,1/EA,6/CV
HC-POPVR-E2-C2	C CounterDisp,3BK,3BL,3PR,3RD

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Handle Height	30.5" - 37.5"
Weight Capacity	350 lb