

# *the Cuff*<sup>®</sup>

*adjustable wrist & ankle weights*

FABRICATION ENTERPRISES INC.  
 TEL: 1-914-345-9300 / 800-431-2830  
 FAX: 1-914-345-9800 / 800-634-5370  
 info@Fab-Ent.com

The Cuff Adjustable Wrist & Ankle Weights



Adult 20 lb Ankle Weight

increase weight as progress improves



adjust the overall weight by changing the number of weight inserts



Pediatric 2 lb Wrist Weight

pediatric and adult sizes available



Adult 4 lb Wrist Weight

- 10-3340-1 Adjustable Wrist Weight, 2 lbs (10 x 0.2 lb inserts)
- 10-3341-1 Adjustable Wrist Weight, 4 lbs (20 x 0.2 lb inserts)
- 10-3331-1 Adjustable Ankle Weight, 5 lbs (10 x 0.5 lb inserts)
- 10-3332-1 Adjustable Ankle Weight, 10 lbs (20 x 0.5 lb inserts)
- 10-3333-1 Adjustable Ankle Weight, 20 lbs (20 x 1.0 lb inserts)
- 10-3345-1 Adjustable Pediatric Wrist Weight, 2 lbs (12 x 0.17 lb inserts)
- 10-3335-1 Adjustable Pediatric Ankle Weight, 2 lbs (12 x 0.17 lb inserts)

Call us for special pricing on weight pairs!



FABRICATION ENTERPRISES INC.  
 PO BOX 1500  
 WHITE PLAINS, NEW YORK 10602 (USA)



AUTHORIZED CE REPRESENTATIVE: RMS UK LTD.  
 28 TRINITY ROAD  
 NAILSEA, SOMESSET BS48 4NU (UK)