


# BASELINE<sup>®</sup> EVALUATION INSTRUMENTS

## Pneumatic Bulb Dynamometer without reset



12-0290



2" squeeze bulb dynamometer measures PSI.

**12-0290:** 0-30 PSI, without max pointer, 0.5 PSI gradations

**12-0292:** 0-15 PSI, without max pointer, 0.5 PSI gradations

This product contains dry natural rubber.

 Fabrication Enterprises Inc.  
250 Clearbrook Rd, Suite 240  
Elmsford, NY 10523 (USA)  
tel: +1-914-345-9300 • 800-431-2830  
fax: +1-914-345-9800 • 800-634-5370  
FabEnt.com

 EC REP  
 AJW Technology Consulting GmbH  
Breite Strasse 3  
40213 Düsseldorf (Germany)