


BASELINE[®] EVALUATION INSTRUMENTS

plastic 6-inch 360° goniometer

REF

12-1002



- permits observation of joint's axis of motion, and range of motion
- 360° head has three scales for ease of use
- primary scale reads in 1° increments
- goniometer arm has a linear scale in inches and centimeters
- goniometer arm remains at measurement until reset

other goniometer styles and sizes available:

12 inch 360°	12-1000*
8 inch 360°	12-1001*
6 inch 360°	12-1002*
6 inch 180° pocket.....	12-1005*
6 inch 360° rulongmeter	12-1006*
3 inch -30° to 120° finger.....	12-1012*

* available standard or with white HiRes[®] strip behind numbers (add HR to end of product code for HiRes[®] model)



FEI
FABRICATION
ENTERPRISES INC.

Fabrication Enterprises Inc.
250 Clearbrook Rd, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300
fax: +1-914-345-9800
FabEnt.com

EC REP



Authorized CE representative:
AJW Technology Consulting GmbH
Königsallee 106
40215 Düsseldorf (Germany)