

**plastic 6-inch 180° pocket
goniometer**

REF

12-1005



- permits observation of joint's axis of motion, and range of motion
- 180° head has three scales for ease of use
- primary scale reads in 1° increments
- goniometer arm has a linear scale in inches and centimeters
- goniometer arm remains at measurement until reset

**other goniometer styles and
sizes available:**

12" 360°	12-1000*
8" 360°	12-1001*
6" 360°	12-1002*
6" 180° pocket	12-1005*
6" 360° rulongmeter	...	12-1006*
3" -30° to 120° finger	..	12-1012*

* available standard or with white HiRes™ strip behind numbers (add HR to end of product code for HiRes™ model)



Fabrication Enterprises Inc.
250 Clearbrook Rd, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
FabEnt.com

EC **REP**



AJW Technology Consulting GmbH
Königsallee 106
40215 Düsseldorf (Germany)