



Overbed Tables

Model 121A



Features:

- High pressure laminate top
- Flush mounted vinyl t-mold edge
- Heavy duty, height adjustable column with a slim and maneuverable base
- Base: Spring assisted height mechanism which locks into place in 1" increments between 29.5" and 44.5"
- Slight upward pressure allows the table to be adjusted without requiring the handle to be activated
- H style base design
- 24" opening
- Standard Chrome plated steel column & base
- Standard Gunstock Walnut Finish
- Optional dual mirror vanity

Specifications:

- Top: 15" x 30"
- Height: 29.5" - 44.5"
- 1.5" Casters
- Weight capacity: 100 lbs evenly distributed
- Warranty: 1 Year
- **123A same table as shown with vanity**

Finish Options



Bannister
Oak

Clear Maple

English Oak

Figured
Mahogany

Fusion
Maple



Gunstock
Walnut

Nepal Teak

Quartered
Oak

Solar Oak

Wild Cherry

Standard Edge Choices



Brown

Opal

Black

Standard Base Finish Choices



Brown

Gray

Opal

Black

Chrome

*Additional costs may apply to finish and trim upgrades