


2650KL - Digital Wheelchair Dual Ramp Scale

The 2650KL offers a remote display head, a medium-sized platform, and folds for convenient storage. The dual ramps make it easy to wheel patients on and off the platform.




2700KL - Digital Wheelchair Dual Ramp Scale with Extra Large Platform

Featuring the largest platform in the wheelchair scale collection, the 2700KL can easily accommodate even the widest wheelchair. It is a good choice in facilities where the scale will remain in one location.

Advantages

- ✓ **EMR Connectivity**
Scale interfaces with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- ✓ **Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's BMI
- ✓ **Versatile Platform**
Wide and low profile platform accommodates wheelchairs and also provides a comfortable weighing experience for bariatric patients
- ✓ **Antimicrobial Keypad**
The keypad features antimicrobial protection
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option, Parent/Child Function) included at no additional cost
- ✓ **Dual Ramps and Remote Display**
Non-slip ramps provide easy on and off access to weighing platform and remote display provides versatility with the ability to be mounted on a wall or placed on a table
- ✓ **Everlock® (patent pending)**
Option that allows users to permanently lock the scale in either LB or KG to improve safety



2650KL

2700KL



Specifications

- Capacity: 1,000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform 2650KL/2650kg (w x d x h):
31 1/2" x 31 1/2" x 2" / 800 mm x 800 mm x 51 mm
- Platform 2700KL/2700KG (w x d x h):
43" x 42" x 2 1/2" / 1092 mm x 1067 mm x 64 mm
- Ramp Size 2650KL/2650KG (w x d): 10" x 31 1/4" / 254 mm x 794 mm
- Ramp Size 2700KL/2700KG (w x d): 12 3/4" x 41 3/4" / 324 mm x 1060 mm
- Product Footprint 2650KL/2650KG (w x h x d):
51 3/4" x 33 5/8" x 5 1/4" / 1314 mm x 854 mm x 133 mm
- Product Footprint 2700KL/2700KG (w x h x d):
69 1/8" x 51 1/2" x 3 3/8" / 1756 mm x 1308 mm x 86 mm
- Product Weight 2650KL/2650KG: 82 lb / 37 kg
- Product Weight 2700KL/2700KG: 187 lb / 85 kg
- Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen
- Connectivity: USB, optional Pelstar[®] wireless technology
- Power Source: 100-240V adapter included, 6 D batteries (not included)
- Functions: LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option, Parent / Child Function, Everlock[®]
- 2 wheels for ease of movement
- 2 year limited warranty

Model	Description
2650KL	Digital wheelchair dual ramp scale, pounds & kilograms
2650KL-BT	2650KL with Pelstar [®] wireless technology
2650KG	Digital wheelchair dual ramp scale, kilograms only
2650KG-BT	2650KG with Pelstar [®] wireless technology
2700KL	Digital wheelchair dual ramp, XL platform, pounds & kilograms
2700KL-BT	2700KL with Pestar [®] wireless technology
2700KG	Digital wheelchair dual ramp, XL platform, kilograms only
2700KG-BT	2700KG with Pelstar [®] wireless technology
SS-2650KL	ScaleSurance 2 year extended warranty for 2650KL, 2650KL-BT, 2650KG & 2650KG-BT*
SS-2700KL	ScaleSurance 2 year extended warranty for 2700KL, 2700KL-BT, 2700KG & 2700KG-BT*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer