



ABD Roll 8" x 20 yd

3001

12/CS

Product Features

Highly absorbent multi-layer dressing

Soft non-woven outer facing wicks moisture into the cellulose core that quickly absorbs and disperses fluid

Moisture resistant barrier prevents fluid strike-through

Sides are sealed to prevent linting and leakage

PreCut

Attributes

Brand Name	Dukal
Size	8" x 20 yd
Sterile	No
Weight	35.64 lbs
Everyday Care	No
Premium Care	No

Packaging

Packaging	12/CS
Case Packaging	12/CS
Total Case Quantity	12
Case Packaging Barcode	(01)20665973005403