

3M™ Tegaderm™ High Performance Foam Adhesive Dressing

90616

large oval

Size: 19 cm x 22.2 cm

Foam: 14 cm x 17.1 cm

90613

medium oval

Size: 14.3 cm x 15.6 cm

Foam: 10 cm x 11 cm

90612

square

Size: 14.3 cm x 14.3 cm

Foam: 10 cm x 10 cm

90619

heel/elbow design

Size: 13.97 cm x 13.97 cm

Foam: 7.62 cm x 7.62 cm

90611

small oval

Size: 10 cm x 11 cm

Foam: 6 cm x 7.6 cm

90614

mini oval

Size: 6.9 cm x 7.6 cm

Foam: 3.1 cm x 3.8 cm

90610

square

Size: 8.8 cm x 8.8 cm

Foam: 5 cm x 5 cm

90615

mini wrap

Size: 6.9 cm x 6.9 cm

Foam: 2.5 cm x 2.5 cm



The mini wrap dressing is constructed of a conformable, absorbent, polyurethane foam pad with a highly breathable, non-waterproof, film backing reinforced with soft cloth tape.

NOTE: Dressings shown are actual size

3M Health Care
3M New Zealand Limited
94 Apollo Drive
Rosedale,
Auckland 0632
Customer Service: 0800 80 81 82
www.3M.com/SkinWoundCare

3M Australia Pty Limited
ABN 90 000 100 096
Building A
1 Rivett Road
North Ryde NSW 2113
1300 363 878
www.3M.com.au/healthcare

© 3M 2012. All rights reserved.
Please recycle. Printed in New Zealand.



Foam Adhesive

With Soft Cloth