

## Convertible Staircase

### MODEL: 4-5080

- Uses steel fasteners and metal braces for extra support and durability
- Pre-assembled handrails and uprights
- One piece handle/rail support system
- Comfort grip, moulded, solid hardwood handrails
- All hardwood plywood base with an ultra-durable UV-cured, topcoat finish
- Easy-grip moulded handrail
- Rigid, bolt-through mounted handrail
- Two, full width, safety treads on each step
- Tough PVC edge tape
- Powder-coated, 2" x 2" steel uprights
- Brace reinforced underside
- Converts from straight to corner arrangement
- Available in two widths
- Requires a minimum 9' ceiling
- 500 lbs. bariatric load capacity under normal use



**IMPORTANT NOTE:** Due to the size of the product, it **MUST** ship to a truck height dock.



California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including vinyl chloride, which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

3

#### MODEL 4-5080-30 SPECIFICATIONS

Platform	30" W x 30" D x 24" H
3 Steps	30" W x 10" D x 6" H
5 Steps	30" W x 10" D x 4" H
Overall Size	62" L

#### MODEL 4-5080-360 SPECIFICATIONS

Platform	36" W x 36" D x 24" H
3 Steps	36" W x 10" D x 6" H
5 Steps	36" W x 10" D x 4" H
Overall Size	62" L