


One-Sided Staircase

MODEL: 4-5501

- Uses steel fasteners and metal braces for support and durability
- Pre-assembled handrails and uprights
- One piece handle/rail support system
- Comfort grip, moulded, solid hardwood handrails
- All hardwood plywood base with an ultra-durable UV-cured, topcoat finish
- Rigid, bolt-through mounted handrail
- Two, full width, safety treads on each step
- Tough PVC edge tape
- Powder-coated, 2" x 2" steel uprights
- Brace reinforced underside
- Small corner staircase for limited space
- Available in 2 widths
- Requires a minimum 9' ceiling
- 500 lbs. bariatric load capacity under normal use



California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including vinyl chloride, which is known to the State of California to cause cancer.
For more information go to www.P65Warnings.ca.gov.

3



MODEL 4-5501-30 SPECIFICATIONS

Platform	30" W x 24" D x 24" H
Each Step	30" W x 10" D x 6" H
Overall Size	34" W x 56" D x 59.5" H

MODEL 4-5501-36 SPECIFICATIONS

Platform	36" W x 24" D x 24" H
Each Step	36" W x 10" D x 6" H
Overall Size	40" W x 56" D x 59.5" H