

## Adjustable One-Sided Staircase

### MODEL: 4-6501

- 1.5 diameter stainless steel handrails with end bumpers
- Fast easy height adjustments with self indexing spring action plungers and numbered height indicator strips (no holes to line up or count)
- Second knob removes any side-to-side play for a rigid upright
- Hardwood plywood base with durable, UV satin, topcoat finish
- Assembled with steel fasteners (No glued wood dowels)
- All edges are banded with tough, Slate Gray, PVC edge material
- 2 full width safety treads on each step
- Available in 30" & 36" widths
- 450 lbs. weight limit (both models)
- Adjusts in 1.5" increments
- Requires a minimum 9' ceiling



California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1



#### MODEL 4-6501-30 SPECIFICATIONS

Platform	30" W x 24" D x 24" H
Each Step	30" W x 10" D x 6" H
Overall Size	34" W x 54" D x 49"- 61.79" H

#### MODEL 4-6501-36 SPECIFICATIONS

Platform	36" W x 24" D x 24" H
Each Step	36" W x 10" D x 6" H
Overall Size	40" W x 54" D x 49"- 61.79" H