

Vinyl Powder Free Examination Glove

This glove is made from the highest quality polyvinyl chloride material and formulation available. Our unique Comfort Formulation allows the glove to fit more comfortably to your hand than traditional vinyl gloves. The 412 Series is an excellent alternative to latex gloves and all of our gloves are produced with DOTP, a safe alternative to the harmful plasticizers found in other vinyl gloves. All of this adds up to a superior vinyl glove in quality, comfort and performance:



- California Proposition 65 Compliant
- Ambidextrous, Beaded Cuff, Single Use Only
- Smooth Finish for Improved Feel / Tactile Sensitivity
- 100% Synthetic Vinyl
- Not Made With Natural Rubber Latex

Tensile Strength - Aged/Unaged

min. 11 MPa / min. 11 MPa

Elongation - Aged/Unaged

min. 350% / min. 350%

Dimensions	Finger	Palm	Cuff	Length	Width
Small	min. 0.08 mm	min. 0.08 mm	min. 0.05 mm	min. 240 mm	min. 85 mm
Medium	min. 0.08 mm	min. 0.08 mm	min. 0.05 mm	min. 240 mm	min. 95 mm
Large	min. 0.08 mm	min. 0.08 mm	min. 0.05 mm	min. 240 mm	min. 105 mm
Extra Large	min. 0.08 mm	min. 0.08 mm	min. 0.05 mm	min. 240 mm	min. 115 mm

Item Numbers	Small	Medium	Large	X-Large
412 Series	412-01	412-02	412-03	412-04

Packaging	Boxes	Cases
packaged by weight	100 Gloves	10 Boxes

