

## H Series, E-Z-Clean, Blood Drawing Chair with ClintonClean™ Flip Arm

### MODEL 6330-P

- E-Z-Clean, one-piece, plastic seat
- Seat is slanted back 6° for comfort and stability
- Height and depth adjustable, ClintonClean™ flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- ClintonClean™ elbow rest
- Height adjustable armrest
- 1-1/4" square, heavy duty, all-welded, tubular steel, neutral gray powder-coated frame (3.18 cm)
- Polished chrome base with black bumper
- Large foot pedal
- 300 lbs. load capacity under normal use (136 kg)

**NOTE:** If padded upholstered arms are selected, use color chart to choose color.

**IMPORTANT NOTE:** Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

**ORDERING NOTE:** MUST add arm suffix number to model number when ordering with 600 arm option, otherwise the chair will be shipped as shown.

**All blood drawing products are non-returnable unless they are unused in unopened, original packaging.**



California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



### MODEL 6330-P SPECIFICATIONS

Seat	18" W x 17" D
Overall	30.75" W x 30" D
Height Range	20.25" — 27.5"
Arm Adjustment	7"

### 600 OPTION



*Padded, Armrests and Flip Arm*

### Standard Upholstery Colors For Optional Arms

