

Lab X Series, Blood Drawing Chair with Padded Flip Arm and Drawer

MODEL 66022-P

- E-Z-Clean, one-piece, black, plastic seat
- Height and depth adjustable, ClintonClean™ flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- ClintonClean™ elbow rest
- ClintonClean™ drawer top
- Height-adjustable, one-piece, plastic drawer
- Extra height reduces bending for attendant
- 1-1/4" square, heavy duty, all-welded, tubular steel, frame (3.18 cm)
- Neutral Gray powder-coated frame
- 300 lbs. load capacity under normal use (136 kg)

NOTE: If padded upholstered arms are selected, use color chart to choose color.

IMPORTANT NOTE: Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

ORDERING NOTE: MUST add arm suffix number to model number when ordering with 600 arm option, otherwise the chair will be shipped as shown.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



MODEL 66022-P SPECIFICATIONS

Seat	18" W x 17" D x 26" H
Overall	40" W x 28" D x 42"H
Arm Height Range	33" — 39"
Arm Adjustment	6"



California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

OPTIONS



600 Padded, Armrest and Flip Arm

Standard Upholstery Colors for Optional Arms

