

67002 Store & Go Phlebotomy Cart

MODEL 67002

- Easy-clean, moulded plastic, rimmed top with convenient plastic storage drawer
- Easy-grip, ClintonClean™ solid plastic handle
- Two, mid-height, storage compartments featuring: one, storage tray with two, plastic, storage bins and a drop-in, large, plastic bin
- Large, stable, 5-leg, chrome base
- Pneumatic, height adjustable cylinder
- White, powder-coated, steel, glove box holder
- Universal attachment rail
- 4" dual-wheel, hooded casters
- Fast, easy assembly with basic tools

NOTE: *Items pictured on top not included*

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.

SPECIFICATIONS

Top: 20.75" x 20" / Including Handle: 24" x 20"

Small Bin Holders: 7.5" x 5" on each side of the post

Small Plastic Bins (2): 6.75" x 3.5"

Medium Bin Holder 10" x 7.25"

Medium Plastic Bin: 9.75" x 7.125"


Plastic Drawer: 17.25" x 7"

Base: 5-leg, 26.5" diameter

Height Range: 31.75" — 42"



California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1