

## 67021 Store & Go Phlebotomy Cart

### MODEL 67021

- ClintonClean™ solid plastic top with adjustable tray brackets
- Easy-grip, ClintonClean™ solid plastic handle
- Two, mid-height storage compartments featuring: one, storage tray with two, plastic, storage bins and a drop-in, large, plastic bin
- Stable, 5-leg, chrome base
- Knob, height adjustable, chrome post
- White, powder-coated steel, glove box holder
- Universal attachment rail
- 3" dual-wheel, hooded casters
- Fast, easy assembly with basic tools

**NOTE:** *Items pictured on top not included*

**All blood drawing products are non-returnable unless they are unused in unopened, original packaging.**

#### SPECIFICATIONS

Top: 20.5" x 16"

Top brackets can accommodate trays from 16" x 15" up to 17.75" x 12"

Small Bin Holders: 7.5" x 5" on each side of the post

Small Plastic Bins (2): 6.75" x 3.5"

Medium Bin holder 10" x 7.25"


Medium plastic bin: 9.75" x 7.125"

Base: 5-leg, 22.5" diameter

Height Range: 26.5" — 37.25"



**California residents please be advised of the following, pursuant to Proposition 65:**

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1