

67022 Store & Go Phlebotomy Cart

MODEL 67022

- Easy-clean, moulded plastic, rimmed top with convenient plastic storage drawer
- Easy-grip, ClintonClean™ solid plastic handle
- Two, mid-height storage compartments featuring: one, storage tray with two, plastic, storage bins, a large, tray section with a large, plastic bin and a lockable* storage bin
- Large, stable, 5-leg, chrome base
- Pneumatic, height adjustable cylinder
- White, powder-coated steel, glove box holder
- Universal attachment rail
- 4" dual-wheel, hooded casters
- Fast, easy assembly with basic tools

NOTE: *Items pictured on top not included*

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.

SPECIFICATIONS

Top: 20.75" x 20" / Including Handle: 24" x 20"

Small Bin Holders: 8.75" x 4" on each side of the post

Small Plastic Bins (2): 6.75" x 3.5"

Large Bin Tray 18.75" x 10.75"

Large Plastic Bin: 15.5" x 9.75"

Lockable Bin: 15.5" x 9.75"

Plastic Drawer: 17.25" x 7"

Base: 5-leg, 26.5" diameter

Height Range: 31.75" — 42"



California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1