

One-Bin Phlebotomy Cart

MODEL 67100

- Gray laminate cabinet construction
- Easy-clean, vacuum-formed, gray, ABS, plastic top with textured finish and raised rim on all sides
- Removable plastic storage bin on easy slide track
16-1/2" L x 10-3/4" W x 5" D (42cm L x 27cm W x 13cm D)
- Easy grip, ClintonClean™ solid plastic handle
- Pneumatic height adjustment
- Stable, 5-leg, aluminum base
- Smooth rolling, 3" ball bearing casters (7.62 cm)
- Optional glove box holder available. Add product
GW-2000 Single White Steel Glove Box Holder

NOTE: Items pictured on top not included

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



MODEL 67100 SPECIFICATIONS

Top	20" W x 20.75" D
Base Diameter	24"
Height Range	33" — 43"
Cabinet	13.5" W x 17.25" D x 7.75"


OPTIONS



GW-2000 Single White Steel Glove Box Holder



Prop 65 Warning

7
 **WARNING:** This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: www.P65Warning.ca.gov.

Laminate Color

Gray