

## Recliner Series Hi-Lo Blood Drawing Chair

### MODEL 6810

- Smooth, fully electric reclining and height adjustment
- Self leveling arms and synchronized footrest that move with chair's seat and backrest
- Keeps patients comfortable and may better accommodate patients who may faint during a procedure
- Hand controller
- All-welded, steel, neutral gray powder-coated frame
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.
- Padded seat, backrest, armrests and footrest
- Armrests flip up out of the way for easy entry and exit
- Adjustable foot levelers
- Easy-clean, durable, stain-resistant vinyl
- 375 lbs. load capacity under normal use (170 kg)

**IMPORTANT NOTE:** Clinton Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.



**All blood drawing products are non-returnable unless they are unused in unopened, original packaging.**

### MODEL 6810 SPECIFICATIONS

Seat	20.5" W x 17" D x 22" – 30" H range
Back	20" W x 32" H with 25° – 75° angle
Footrest	19.5" W x 14" D with 5° – 90° angle
Overall	35" W x 50" D x 35" retracted/66" extended
Arm Adjustment	21" between arms



California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

### Standard Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black