



Babinski Percussion Hammer, Adjustable

7018

1/EA

Product Features

Measures nerve response

Chrome plated telescoping handle adjusts for 6 ½" to 15"

Head may be used in both the horizontal or vertical positions

Attributes

Brand Name	Tech-Med
Material	Chrome Plated Steel
Sterile	No
Weight	0.33069 lbs
Everyday Care	No
Premium Care	No

Packaging

Packaging	1/EA
Case Packaging	1/EA
Total Case Quantity	1
Item Packaging Barcode	(01)10665973022779
Case Packaging Barcode	(01)20665973022776