

General, Flat Top, Ultrasound Table



MODEL 80061

- One piece flat top with firm foam padding
- Trendelenburg/Reverse Trendelenburg plus or minus 15 °
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Industry leading 5-year, limited warranty
- Optional side rails & casters available
- 34" wide model available
- Model 80061: 500 lbs load capacity under normal use
- Model 80061-X: 600 lbs load capacity under normal use



Trendelenburg



Reverse Trendelenburg




Optional Side Rails & Casters

80061 SPECIFICATIONS

	length	width	height range
overall	72"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "

80061-X Extra-Wide 34" width model

California residents please be advised of the following, pursuant to Proposition 65:

1
 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

OPTIONS

Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set

080	length	width
	92"	2"
	233.67 cm	5 cm



Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

098



Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087



Paper Cutter

- Clear plastic with elastic strap and buckle

040

