

## Cabinet Style, Treatment Table with Drawers & Doors



### MODEL 9105-30

- Features 3 drawers, adjustable backrest and 1 storage compartment with 1 adjustable shelf
- Easy-clean laminate surfaces
- Flat panel drawers with all steel sides
- Adjustable doors with concealed soft-close hinges
- Ships assembled
- Seamless rounded corner top
- Heavyweight, knit-backed upholstery
- 2" firm foam padding (5 cm)
- 425 lbs. load capacity under normal use

### SPECIFICATIONS

9105-30

72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

### Options



### Laminate Colors



### Upholstery Colors

